

HOOK2 Rabbit Polyclonal Antibody



CAB8210

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

100kDa

Calculated MW:

83kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

29911

Uniprot

Q96ED9

Synonyms:

HOOK2; HK2

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 438-717 of human HOOK2 (NP_001093646.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
