

Product Information

Product SKU:	CAB8225	Gene ID:	54910	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	110kDa	Conjugate:	Unconjugated
Calculated MW:	93kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable chemorepellent activity and semaphorin receptor binding activity. Involved in muscle cell differentiation and positive regulation of stress-activated MAPK cascade. Located in extracellular space.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	SEMAF; SEMAI; SEMACL1; M-SEMA-F; SEMA4C
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 749-833 of human SEMA4C (NP_060259.4).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.