

Product Information

Product SKU:	CAB8257	Gene ID:	79132	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	77kDa	Conjugate:	Unconjugated
Calculated MW:	77kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables double-stranded RNA binding activity; single-stranded RNA binding activity; and zinc ion binding activity. Involved in negative regulation of defense response and negative regulation of type I interferon production. Predicted to be active in cytoplasm.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Synonyms:	LGP2; RLR-3; D11LGP2; D11lgp2e; DHX58
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 389-678 of human DHX58 (NP_077024.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.