

RNF219 Rabbit Polyclonal Antibody



CAB8259

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

81kDa

Calculated MW:

81kDa

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

79596

Uniprot

Q5W0B1

Synonyms:

RNF219; C13orf7


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 600 to the C-terminus of human RNF219 (NP_078822.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - RNF219 Polyclonal Antibody (CAB8259)