

RSAD2 Rabbit Polyclonal Antibody



CAB8271

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

46KDa

Calculated MW:

42kDa

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

91543

Uniprot

Q8WVG1

Synonyms:

RSAD2; 2510004L01Rik; cig33; cig5; vig1


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 42-361 of human RSAD2 (NP_542388.2).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using RSAD2 antibody (CAB8271) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (CABM00021). Exposure time: 180s.



Immunofluorescence analysis of U2OS cells using RSAD2 antibody (CAB8271) at dilution of 1:100. Blue: DAPI for nuclear staining.