

CAB8280

Product Information

Product SKU:	CAB8280	Gene ID:	147906	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	65kDa	Conjugate:	Unconjugated
Calculated MW:	65kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable delta-catenin binding activity; protein kinase A binding activity; and protein kinase C binding activity. Involved in negative regulation of canonical Wnt signaling pathway and negative regulation of cell growth. Predicted to be active in cytoplasm.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: RRR1; DAPPER3; DACT3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human DACT3 (NP_659493.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.