

Product Information

| | | | | | |
|--------------|---------|---------------|--------|-------------|-------------|
| Product SKU: | CAB8287 | Gene ID: | 252983 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse |

Additional Information

| | | | |
|----------------|-------|------------|--------------|
| Observed MW: | 70kDa | Conjugate: | Unconjugated |
| Calculated MW: | 62kDa | Isotype: | IgG |

Immunogen Information

Background: Enables syntaxin binding activity. Involved in several processes, including positive regulation of cell cycle G1/S phase transition; positive regulation of keratinocyte proliferation; and protein stabilization. Located in extracellular exosome.

| | |
|-----------------------|--|
| Recommended Dilution: | WB,1:500 - 1:2000 IHC-P,1:50 - 1:100 |
| Synonyms: | Synip; STXBP4 |
| Purification Method: | Affinity purification |
| Immunogen: | Recombinant fusion protein containing a sequence corresponding to amino acids 344-553 of human STXBP4 (NP_848604.3). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |