

SIK2 Rabbit Polyclonal Antibody



CAB8321

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

129kDa

Calculated MW:

103kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

23235

Uniprot

Q9H0K1

Synonyms:

SIK2; LOH11CR1; QIK; SIK-2; SNF1LK2


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 827-926 of human SIK2 (NP_056006.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using SIK2 antibody (CAB8321) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 30s.