

Product Information

Product SKU:	CAB8376	Gene ID:	148581	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	42kDa	Conjugate:	Unconjugated
Calculated MW:	38kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable ubiquitin conjugating enzyme activity. Predicted to be involved in DNA repair; proteasome-mediated ubiquitin-dependent protein catabolic process; and protein ubiquitination. Predicted to be part of HULC complex.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:100
Synonyms:	UBE2U
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-226 of human UBE2U (NP_689702.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.