

MEX3C Rabbit Polyclonal Antibody



CAB8393

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

73kDa

Calculated MW:

69kDa

Applications:

WB

Reactivity:

Human

Protein Background

This gene encodes a member of a family of proteins with two K homology (KH) RNA-binding domains and a C-terminal RING-finger domain. The protein interacts with mRNA via the KH domains, and the protein shuttles between the nucleus and cytoplasm. Polymorphisms in this gene may contribute to hypertension.

Immunogen information

Gene ID:

51320

Uniprot

Q5U5Q3

Synonyms:

MEX3C; BM-013; MEX-3C; RKHD2; RNF194; Mex3C

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 484-659 of human MEX3C (NP_057710.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
