

# CDO1 Rabbit Polyclonal Antibody



CAB8402

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

23kDa

### Calculated MW:

22kDa

### Applications:

WB

### Reactivity:

Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

1036

#### Uniprot

Q16878

#### Synonyms:

CDO1; CDO-I

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CDO1 (NP\_001792.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---