

Product Information

Product SKU:	CAB8407	Gene ID:	1380	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	150kDa	Conjugate:	Unconjugated
Calculated MW:	113kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a membrane protein, which functions as a receptor for Epstein-Barr virus (EBV) binding on B and T lymphocytes. Genetic variations in this gene are associated with susceptibility to systemic lupus erythematosus type 9 (SLEB9). Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:100 IF/ICC,1:50 - 1:100

Synonyms: CR; C3DR; CD21; CVID7; SLEB9; CR2/CD21

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 300-590 of human CR2/CD21 (NP_001006659.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.