

Product Information

Product SKU:	CAB8423	Gene ID:	2790	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	14kDa	Conjugate:	Unconjugated
Calculated MW:	7kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable G-protein beta-subunit binding activity. Predicted to be involved in G protein-coupled receptor signaling pathway. Predicted to be located in plasma membrane. Predicted to be part of heterotrimeric G-protein complex.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	GNG10
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-68 of human GNG10 (NP_001017998.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.