

Product Information

Product SKU:	CAB8502	Gene ID:	56981	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	58kDa	Conjugate:	Unconjugated
Calculated MW:	58kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable chromatin binding activity. Involved in several processes, including negative regulation of cell growth; positive regulation of fibroblast apoptotic process; and regulation of transcription, DNA-templated. Located in cytosol and nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	PFM8; PRDM11
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 228-477 of human PRDM11 (NP_001243625.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.