

Product Information

Product SKU:	CAB8511	Gene ID:	79736	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	50kDa	Conjugate:	Unconjugated
Calculated MW:	42kDa	Isotype:	IgG

Immunogen Information

Background: Enables DNA polymerase processivity factor activity. Involved in mitochondrial transcription and oxidative phosphorylation. Located in mitochondrial nucleoid. Part of ribonucleoprotein complex.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200
Synonyms:	C17orf42; TEFM
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 36-210 of human TEFM (NP_078959.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.