

Product Information

Product SKU:	CAB8576	Gene ID:	23130	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	213kDa	Conjugate:	Unconjugated
Calculated MW:	213kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable phosphatidylinositol-3-phosphate binding activity. Involved in autophagosome assembly. Predicted to be located in endoplasmic reticulum membrane; lipid droplet; and phagophore assembly site membrane. Predicted to be active in phagophore assembly site. Predicted to be extrinsic component of membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	BLTP4A; ATG2A
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 700-800 of human ATG2A (NP_055919.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.