

Product Information

Product SKU:	CAB8594	Gene ID:	114971	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	23kDa	Conjugate:	Unconjugated
Calculated MW:	23kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable phosphatidylglycerophosphatase activity and phosphatidylinositol-4,5-bisphosphate 5-phosphatase activity. Involved in regulation of intrinsic apoptotic signaling pathway. Located in mitochondrion.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	MOSP; PLIP; DUSP23; PNAS-129; PTPMT1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-201 of human PTPMT1 (NP_783859.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.