

PTPMT1 Rabbit Polyclonal Antibody



CAB8594

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

23kDa

Calculated MW:

15kDa/16kDa/22kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

114971

Uniprot

Q8WUK0

Synonyms:

PTPMT1; DUSP23; MOSP; PLIP; PNAS-129

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-201 of human PTPMT1 (NP_783859.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
