

SLC16A11 Rabbit Polyclonal Antibody



CAB8599

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

48kDa

Calculated MW:

48kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

216867

Uniprot

Synonyms:

SLC16A11; MCT11


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse SLC16A11 (NP_694721.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - SLC16A11 Polyclonal Antibody (CAB8599)