

HNRNPCL1 Rabbit Polyclonal Antibody



CAB8606

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

32kDa

Applications:

IHC

Reactivity:

Human, Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

343069

Uniprot

O60812

Synonyms:

HNRNPCL1; HNRPCL1


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human HNRNPCL1 (NP_001013653.1).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded rat kidney using HNRNPCL1 antibody (CAB8606) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human stomach using HNRNPCL1 antibody (CAB8606) at dilution of 1:100 (40x lens).