

Product Information

Product SKU:	CAB8608	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	51kDa	Conjugate:	Unconjugated
Calculated MW:	-	Isotype:	IgG

Immunogen Information

**Background:** Hand, foot and mouth disease is caused by different viruses from the enterovirus family, mainly coxsackievirus A16 and enterovirus 71 (EV71) from the human enterovirus A family. Severe disease and neurological complications are associated more often with EV71 infection, and can lead occasionally to fatal brain stem encephalitis in young children.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	-
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of Enterovirus EV71 3D
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.