

p150 CAF1 Rabbit Monoclonal Antibody
CAB8725



Product Information

Product SKU:	CAB8725	Gene ID:	10036	Size:	20uL, 100uL
Clone No:	ARC2011	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	150kDa	Conjugate:	Unconjugated
Calculated MW:	107kDa	Isotype:	IgG

Immunogen Information

Background: Chromatin assembly factor I (CAF1) is a nuclear complex consisting of p50, p60 (CHAF1B; MIM 601245), and p150 (CHAF1A) subunits that assembles histone octamers onto replicating DNA in vitro (Kaufman et al., 1995 [PubMed 7600578]).

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	CAF1; P150; CAF-1; CAF1B; CAF1P150; p150 CAF1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human p150 CAF1 (NP_005474.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.