

Product Information

Product SKU:	CAB8731	Gene ID:	84951	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	77kDa	Conjugate:	Unconjugated
Calculated MW:	77kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable actin binding activity. Involved in protein localization. Located in focal adhesion.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	CTEN; PP14434; TNS4
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 380-520 of human TNS4 (NP_116254.4).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.