

PMP2 Rabbit Polyclonal Antibody



CAB8851

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

15kDa

Calculated MW:

14kDa

Applications:

WB

Reactivity:

Mouse

Protein Background

Immunogen information

Gene ID:

5375

Uniprot

P02689

Synonyms:

PMP2; FABP8; M-FABP; MP2; P2

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-132 of human PMP2 (NP_002668.1).

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of mouse small intestine, using PMP2 antibody (CAB8851) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 60s.