

# C1QTNF4 Rabbit Polyclonal Antibody



CAB8954

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

38kDa

### Calculated MW:

35kDa

### Applications:

WB

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

114900

#### Uniprot

Q9BXJ3

#### Synonyms:

C1QTNF4; CTRP4; ZACRP4

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 250-329 of human C1QTNF4 (NP\_114115.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---