

Product Information

Product SKU:	CAB9061	Gene ID:	28985	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	20kDa	Conjugate:	Unconjugated
Calculated MW:	21kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable translation initiation factor activity. Involved in IRES-dependent viral translational initiation; formation of translation preinitiation complex; and ribosome disassembly. Located in cytosol and plasma membrane. Colocalizes with cytosolic small ribosomal subunit.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Synonyms:	MCT1; MCT-1; MCTS1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-181 of human MCTS1 (NP_054779.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.