

# PDCD4 Rabbit Monoclonal Antibody



CAB9068

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

52KDa

### Calculated MW:

52kDa

### Applications:

WB IHC

### Reactivity:

Human, Mouse, Rat

## Protein Background

This gene is a tumor suppressor and encodes a protein that binds to the eukaryotic translation initiation factor 4A1 and inhibits its function by preventing RNA binding. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2010]

## Immunogen information

### Gene ID:

27250

### Uniprot

Q53EL6

### Synonyms:

H731

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50  
- 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:


A synthesized peptide derived from human PDCD4

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

## Product Images

---



Western blot - PDCD4 Rabbit mAb (CAB9068)