

Product Information

Product SKU:	CAB9143	Gene ID:	8662	Size:	20uL, 100uL
Clone No:	ARC1448	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	116kDa	Conjugate:	Unconjugated
Calculated MW:	93kDa	Isotype:	IgG

Immunogen Information

Background:	Enables RNA binding activity. Contributes to translation initiation factor activity. Involved in several processes, including IRES-dependent viral translational initiation; translational initiation; and viral translational termination-reinitiation. Located in extracellular exosome. Part of eukaryotic translation initiation factor 3 complex.
Recommended Dilution:	WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200
Synonyms:	PRT1; EIF3S9; EIF3-ETA; EIF3-P110; EIF3-P116; eIF3B
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 101-200 of human eIF3B (NP_003742.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.