

CAB9161

Product Information

Product SKU:	CAB9161	Gene ID:	23279	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	160kDa	Conjugate:	Unconjugated
Calculated MW:	162kDa	Isotype:	IgG

Immunogen Information

Background: A structural constituent of nuclear pore. Involved in mRNA export from nucleus and nephron development. Part of nuclear pore outer ring. Colocalizes with kinetochore. Implicated in nephrotic syndrome type 19.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	NPHS19; NUP160
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 40-150 of human NUP160 (NP_056046.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.