

Product Information

Product SKU:	CAB9193	Gene ID:	6829	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	160kDa	Conjugate:	Unconjugated
Calculated MW:	121kDa	Isotype:	IgG

Immunogen Information

Background: Enables enzyme binding activity and protein heterodimerization activity. Involved in positive regulation of macroautophagy; regulation of RNA metabolic process; and transcription elongation from RNA polymerase II promoter. Located in nucleoplasm. Part of DSIF complex.

Recommended Dilution:	WB,1:200 - 1:2000 IP,0.5µg-4µg antibody for 200µg-400µg extracts of whole cells
Synonyms:	SPT5; SPT5H; Tat-CT1; SUPT5H/SPT5
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 125-425 of human SUPT5H/SPT5 (NP_003160.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.