

ALKBH1 Rabbit Monoclonal Antibody



CAB9221

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

43KDa

Calculated MW:

43kDa

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Protein Background

This gene encodes a homolog to the E. coli alkB gene product. The E. coli alkB protein is part of the adaptive response mechanism of DNA alkylation damage repair. It is involved in damage reversal by oxidative demethylation of 1-methyladenine and 3-methylcytosine. [provided by RefSeq, Jul 2008]

Immunogen information

Gene ID:

8846

Uniprot

Q13686

Synonyms:

ABH; ABH1; ALKBH; alkB; hABH

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

A synthesized peptide derived from human ALKBH1

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images



Western blot - ALKBH1 Rabbit mAb (CAB9221)