

CAB9232

Product Information

Product SKU:	CAB9232	Gene ID:	55773	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	78kDa	Conjugate:	Unconjugated
Calculated MW:	78kDa	Isotype:	IgG

Immunogen Information

Background: Involved in brain development; retrograde transport, endosome to Golgi; and vesicle tethering to Golgi. Located in cytoplasmic vesicle and trans-Golgi network. Colocalizes with WASH complex. Implicated in pontocerebellar hypoplasia.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: PCH11; NS4ATP1; TBC1D23

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 250-600 of human TBC1D23 (NP_001186127.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.