

# TEKT4 Rabbit Polyclonal Antibody



CAB9241

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

50kDa

### Applications:

WB

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:1000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

150483

#### Uniprot

Q8WW24

#### Synonyms:

TEKT4; tektin-4

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human TEKT4 (NP\_653306.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---