

USP22 Rabbit Monoclonal Antibody

CAB9261



Product Information

Product SKU:	CAB9261	Gene ID:	23326	Size:	20uL, 100uL
Clone No:	ARC1498	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	60kDa	Conjugate:	Unconjugated
Calculated MW:	60kDa	Isotype:	IgG

Immunogen Information

Background:	Enables enzyme binding activity; nuclear receptor coactivator activity; and thiol-dependent deubiquitinase. Contributes to H4 histone acetyltransferase activity. Involved in positive regulation of mitotic cell cycle; positive regulation of transcription, DNA-templated; and protein modification by small protein conjugation or removal. Part of SAGA complex.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	USP3L; USP22
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-175 of human USP22 (Q9UPT9).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.