

CAB9263

Product Information

Product SKU:	CAB9263	Gene ID:	1175	Size:	20uL, 100uL
Clone No:	ARC1500	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	17kDa	Conjugate:	Unconjugated
Calculated MW:	17kDa	Isotype:	IgG

Immunogen Information

Background: One of two major clathrin-associated adaptor complexes, AP-2, is a heterotetramer which is associated with the plasma membrane. This complex is composed of two large chains, a medium chain, and a small chain. This gene encodes the small chain of this complex. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: AP17; FBH3; HHC3; FBHOK; CLAPS2; AP2S1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human AP2S1 (P53680).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.