

CAB9476

Product Information

Product SKU:	CAB9476	Gene ID:	7547	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	51kDa	Conjugate:	Unconjugated
Calculated MW:	51kDa	Isotype:	IgG

Immunogen Information

**Background:** This gene encodes a member of the ZIC family of C2H2-type zinc finger proteins. This nuclear protein probably functions as a transcription factor in early stages of left-right body axis formation. Mutations in this gene cause X-linked visceral heterotaxy, which includes congenital heart disease and left-right axis defects in organs.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	HTX; HTX1; ZNF203; VACTERLX; ZIC3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human ZIC3 (NP_003404.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.