

TBC1D1 Rabbit Polyclonal Antibody



CAB9501

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

133kDa

Calculated MW:

131kDa/133kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

23216

Uniprot

Q86TI0

Synonyms:

TBC1D1; TBC; TBC1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 520-710 of human TBC1D1 (NP_055988.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
