

Product Information

Product SKU:	CAB9538	Gene ID:	10598	Size:	20uL, 100uL
Clone No:	ARC1623	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	42kDa	Conjugate:	Unconjugated
Calculated MW:	38kDa	Isotype:	IgG

Immunogen Information

Background: Enables ATPase activator activity; Hsp90 protein binding activity; and chaperone binding activity. Involved in positive regulation of ATPase activity. Located in cytosol.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	p38; AHA1; hAha1; C14orf3
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 239-338 of human AHA1 (O95433).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.