

---

## Product Information

<b>Product SKU:</b>	CAB9555	<b>Gene ID:</b>	84083	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	123kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	123kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** Enables ATP-dependent DNA/DNA annealing activity; K63-linked polyubiquitin modification-dependent protein binding activity; and endodeoxyribonuclease activity. Involved in several processes, including DNA metabolic process; DNA rewinding; and negative regulation of DNA recombination. Located in nuclear replication fork and nucleoplasm.

**Recommended Dilution:** WB,1:1000 - 1:3000

**Synonyms:** AH2; 4933425L19Rik; ZRANB3

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 818-1077 of human ZRANB3 (NP\_001273497.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.