

Product Information

Product SKU:	CAB9582	Gene ID:	10713	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	65kDa	Conjugate:	Unconjugated
Calculated MW:	65kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable thiol-dependent deubiquitinase and zinc ion binding activity. Involved in spliceosomal complex assembly. Located in nucleoplasm. Part of U4/U6 x U5 tri-snRNP complex.

Recommended Dilution:	WB,1:1000 - 1:2000
Synonyms:	65K; SAD1; CGI-21; HSPC332; SNRNP65; USP39
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 336-565 of human USP39 (NP_006581.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.