

Product Information

Product SKU:	CAB9585	Gene ID:	11345	Size:	20uL, 100uL
Clone No:	ARC1648	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	16kDa	Conjugate:	Unconjugated
Calculated MW:	14kDa	Isotype:	IgG

Immunogen Information

Background: Enables ubiquitin protein ligase binding activity. Involved in negative regulation of proteasomal protein catabolic process and protein localization to endoplasmic reticulum. Located in Golgi membrane and autophagosome membrane.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	ATG8; GEF2; ATG8C; GEF-2; GATE16; GATE-16; FLC3A; GABARAPL2
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 38-117 of human GABARAPL2 (P60520).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.