

CAB9617

---

## Product Information

<b>Product SKU:</b>	CAB9617	<b>Gene ID:</b>	27341	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	37kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	32kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** Enables RNA binding activity. Predicted to be involved in rRNA processing and ribosomal small subunit assembly. Predicted to act upstream of or within blastocyst formation. Predicted to be located in nucleoplasm. Predicted to be part of CURI complex and UTP-C complex. Implicated in primary autosomal recessive microcephaly.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** Rrp7; CGI-96; MCPH28; BK126B4.3; RRP7A

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human RRP7A (NP\_056518.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.