

Cytokeratin 9 (KRT9) Rabbit Monoclonal Antibody

CAB9621



Product Information

Product SKU:	CAB9621	Gene ID:	3857	Size:	20uL, 100uL
Clone No:	ARC1664	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	62kDa	Conjugate:	Unconjugated
Calculated MW:	62kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes the type I keratin 9, an intermediate filament chain expressed only in the terminally differentiated epidermis of palms and soles. Mutations in this gene cause epidermolytic palmoplantar keratoderma.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	K9; CK-9; EPPK; Cytokeratin 9 (KRT9)
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 254-412 of human Cytokeratin 9 (KRT9) (P35527).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.