

CAB9634

Product Information

Product SKU:	CAB9634	Gene ID:	1656	Size:	20uL, 100uL
Clone No:	ARC1673	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	54kDa	Conjugate:	Unconjugated
Calculated MW:	54kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the DEAD box protein family. The protein is an RNA helicase found in P-bodies and stress granules, and functions in translation suppression and mRNA degradation. It is required for microRNA-induced gene silencing. Multiple alternatively spliced variants, encoding the same protein, have been identified.

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

Synonyms: P54; RCK; HLR2; IDDILF; Rck/p54; DDX6

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 101-200 of human DDX6 (P26196).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.