

CAB9644

Product Information

Product SKU:	CAB9644	Gene ID:	5360	Size:	20uL, 100uL
Clone No:	ARC1676	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	80kDa	Conjugate:	Unconjugated
Calculated MW:	55kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is one of at least two lipid transfer proteins found in human plasma. The encoded protein transfers phospholipids from triglyceride-rich lipoproteins to high density lipoprotein (HDL). In addition to regulating the size of HDL particles, this protein may be involved in cholesterol metabolism. At least two transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: BPIFE; HDLCQ9; PLTP

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PLTP (P55058).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.