

CAB9666

Product Information

Product SKU:	CAB9666	Gene ID:	56993	Size:	20uL, 100uL
Clone No:	ARC1688	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	18kDa	Conjugate:	Unconjugated
Calculated MW:	16kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is an integral membrane protein of the mitochondrial outer membrane. The encoded protein interacts with TOMM20 and TOMM40, and forms a complex with several other proteins to import cytosolic preproteins into the mitochondrion.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IP,0.5µg-4µg antibody for 200µg-400µg extracts of whole cells

Synonyms: 1C9-2; TOM22; MST065; MSTP065; TOMM22

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TOMM22 (Q9NS69).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.