

**GSTK1 Rabbit Monoclonal Antibody**

**CAB9667**



**Product Information**

<b>Product SKU:</b>	CAB9667	<b>Gene ID:</b>	373156	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC1689	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

**Additional Information**

<b>Observed MW:</b>	25kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	25kDa	<b>Isotype:</b>	IgG

**Immunogen Information**

**Background:** This gene encodes a member of the kappa class of the glutathione transferase superfamily of enzymes that function in cellular detoxification. The encoded protein is localized to the peroxisome and catalyzes the conjugation of glutathione to a wide range of hydrophobic substates facilitating the removal of these compounds from cells. Alternative splicing results in multiple transcript variants.

<b>Recommended Dilution:</b>	WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200
<b>Synonyms:</b>	GST; GST13; hGSTK1; GSTK1-1; GST13-13; GST 13-13; GSTK1
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human GSTK1 (Q9Y2Q3).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.