

Product Information

Product SKU:	CAB9689	Gene ID:	23016	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	37kDa	Conjugate:	Unconjugated
Calculated MW:	32kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable 3'-5'-exoribonuclease activity and RNA binding activity. Predicted to be involved in RNA metabolic process. Part of exosome (RNase complex).

Recommended Dilution:	WB,1:200 - 1:2000
Synonyms:	p8; EAP1; RRP42; Rrp42p; hRrp42p; EXOSC7
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-291 of human EXOSC7 (NP_055819.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.