

PTPRZ1 Rabbit Polyclonal Antibody



CAB9766

Product Information

Size:

100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

162kDa/253kDa/254kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:1000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a member of the receptor protein tyrosine phosphatase family. Expression of this gene is restricted to the central nervous system (CNS), and it may be involved in the regulation of specific developmental processes in the CNS. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

Immunogen information

Gene ID:

5803

Uniprot

P23471

Synonyms:

PTPRZ1; HPTPZ; HPTPzeta; PTP-ZETA; PTP18; PTPRZ; PTPZ; R-PTP-zeta-2; RPTPB; RPTPbeta

Immunogen:

A synthetic peptide of human PTPRZ1

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
