

# BZR2 Rabbit Polyclonal Antibody



CAB9794

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to Figures

### Calculated MW:

36kDa

### Applications:

WB

### Reactivity:

Arabidopsis thaliana

## Protein Background

### Immunogen information

#### Gene ID:

838518

#### Uniprot

#### Synonyms:

BES1; BZR2; 107 PROTEIN; BRASSINAZOLE-RESISTANT 2; BRI1-EMS-SUPPRESSOR 1; F18O14.7; F18O14\_7

#### Immunogen:

Recombinant protein of Arabidopsis thaliana BZR2

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---